






MDL Supplies Request Form

Orders are completed within 1 – 5 business days. Please send this form to mdl@wichita.edu

MDL relies on ordering facilities to acknowledge and track expiration dates. MDL will reject specimens received in expired tubes or bottles.

<u>Item Name</u>	<u>Picture</u>	<u>Quantity</u>
Aerobic Blood Culture Bottles (FA Plus) – Pack of 10		
Anaerobic Blood Culture Bottles (FN Plus) – Pack of 10		
Biohazard Transport Bags – Pack of 100		
ESwab – Pack of 10		
ESwab with Mini-Tip Flocked Swab – Pack of 10		
Pediatric Blood Culture Bottles (PF Plus) – Pack of 10		
RIC Nasopharyngeal Collection Kit – Kit of 10		

RIC Oropharyngeal Collection Kit – Kit of 10		
SalivaDirect Collection Kit – Kit of 10		
Sterile Container – Pack of 10		
Urine Culture Kit – Kit of 10		

Contact Information

Organization Name:

Requester's Name:








Requester's Email Address:

Additional Notes:

Orders are completed within 1 – 5 business days. Please send this form to mdl@wichita.edu

MDL tracks test volume trends. If testing volume does not match supply requests, MDL will reach out to ordering facilities prior to fulfilling orders.

Collection Supplies Guide

Supply Item & Description	Supply Picture	Associated Tests
Blood Culture Bottles – BioMerieux Aerobic Blood Culture Bottle (FA Plus), Anaerobic Blood Culture Bottle (FN Plus), Pediatric Blood Culture Bottle (PF Plus)		Blood Culture
ESwab – Copan Diagnostics ESwab Collection and Transport System		Genital Culture Group A Screen by Culture Group B Screen by Culture MRSA Screen by Culture <i>Neisseria gonorrhoeae</i> Culture Upper Respiratory Culture Wound Culture
ESwab – Mini-Tip Flocked Swab		Genital Culture (Urogenital Source) <i>Neisseria gonorrhoeae</i> Culture (Urogenital Source)
Saliva Collection Kit – 5 mL Eppendorf Tube		SalivaDirect – SARS-CoV-2 by RT-PCR
Sterile Container – 50 mL Sterile Tube		Body Fluid Culture Cerebral Spinal Fluid Culture Lower Respiratory Culture Tissue Culture
RIC Collection Kit – Viral Transport Media with Flocked Swab or Foam Swab		RIC – RSV A/B, Influenza A/B, SARS-CoV-2 panel by RT-PCR
Urine Collection Kit – BD Vacutainer Plus C&S Preservative Tubes and a Urine Transfer Straw		Urine Culture